

Digitimer DS4

Bi-phasic Stimulus Isolator

OPERATOR'S MANUAL



"Digitimer" is a registered trade mark of Digitimer Limited

Product Registration

Please take time out to register your new product.

Details are included in the Introduction.

You can even do it online at:-

www.digitimer.com/register

Table of Contents

INDEX OF FIGURES & TABLES.....	4
Introduction.....	5
Precautions and Warnings.....	6
Operator's Manual.....	6
Explosion and Fire.....	6
Damage.....	6
Moisture.....	6
Electrical Interference.....	6
Contact Address.....	6
Servicing & Maintenance.....	6
Environmental Considerations.....	6
Product Registration.....	7
Why Register your purchase?.....	7
How to Register your purchase.....	7
Product announcement mailing list.....	7
Unpacking.....	7
Included Items.....	7
Optional Accessories.....	7
Declaration of Conformity.....	8
Functions & Features.....	9
Overview of the Controls/Connections.....	9
Power Connection.....	9
Power Zone.....	9
Gate Zone.....	10
Input Zone.....	10
Output Zone.....	10
Output Characteristics.....	11
Checking Actual Output Current.....	12
Output Capacitance and Risetime.....	12
How Constant is the Constant Current Output?.....	12
Frequency Response.....	13
Inactivity Sensor.....	13
Internal Jumper Selection.....	14
Battery Life.....	15

Battery Testing.....	16
Battery Replacement.....	16
Rack Mounting.....	16
References.....	17
DS4 Specifications.....	18
Warranty.....	19
Obtaining Warranty Service.....	19
Product change or discontinuation.....	19
Operator Notes.....	20

Index of Figures & Tables

Figure 1 The DS4 Bi-phasic Stimulus Isolator.....	5
Figure 2 Positioning of the DC power supply socket on the DS4.....	9
Figure 3 The POWER zone of the DS4 front panel.....	9
Figure 4 The GATE zone of the DS4 front panel.....	10
Figure 5 The INPUT zone of the DS4 front panel.....	10
Figure 6 The OUTPUT zone of the DS4 front panel.....	11
Figure 7 Graph of deliverable current (hatched area) depending upon load.....	11
Figure 8 Method to monitor current actually delivered.....	12
Figure 9 Frequency Response Characteristics.....	13
Figure 10 Effect of the inactivity sensor on stimulus initiation.....	14
Figure 11 Location of the internal jumpers.....	14
Table 1 Expected "Control" battery life.....	15
Table 2 Expected "Stimulus" battery life.....	15
Figure 12 DS4 battery compartment labelling.....	16

Introduction

The DS4 has been developed to meet the needs of laboratory scientists who require a stimulus isolator that can output a bi-phasic isolated constant current stimulus in response to an external analogue command voltage signal, provided by a computer DAC via software or a signal generator. Such a requirement is already met by our NeuroLog System in the form of the NL512 Biphasic Buffer and NL800A Stimulus Isolators, but the DS4 provides our first standalone device to meet this need.

The DS4 accepts a variety of voltage input ranges ($\pm 1V$, $\pm 2.5V$, $\pm 5V$ and $\pm 10V$) and produces a constant current stimulus output in 4 overlapping ranges ($\pm 10\mu A$, $\pm 100\mu A$, $\pm 1mA$ and $\pm 10mA$) from a compliance voltage of approximately $\pm 44V$. In addition, the DS4 has a GATE input which allows multiple DS4's to be connected to a single analogue voltage source, with each DS4 being digitally enabled, separately.

One of the problems with stimulators that make use of an external voltage source to define a stimulus waveform is that small offsets or noisy baseline signals from the DAC's used to drive them can result in unwanted battery drain or perhaps worse, low amplitude stimulation. The DS4 uses a special "inactivity sensor" to monitor the input voltage and disable the DS4 output if this voltage falls below 0.15% of the full scale value (in -ve or +ve directions) for a user selectable time period of 200ms, 500ms, 1s or 2s. Unlike other devices which only produce an output when the input voltage exceeds a threshold value, this "inactivity sensor" reduces battery usage and damaging "leak currents" during infrequent stimulation, while at the same time maintaining low levels of zero crossing distortion for repetitive waveforms.

The DS4 uses an external power supply to power the input control circuitry and readily available/inexpensive batteries to provide the opto-isolated stimulus voltage source. Batteries are only used when a stimulus is commanded by an input voltage waveform.



Figure 1 The DS4 Bi-phasic Stimulus Isolator.

Precautions and Warnings

Operator's Manual

Carefully study this Operators Manual before using the DS4 Bi-phasic Stimulus Isolator.

Explosion and Fire

The DS4 must not be used in an explosive atmosphere or volatile atmosphere.

Damage

The DS4 and/or any accessories must not be used if there are any signs of external damage.

Moisture

The DS4 and/or any accessories must not be used if any parts are wet or damp.

Electrical Interference

This unit has been fully tested for European (CE) EMC conformity. This unit should NOT be used near radio transmitters. If any 'strange' behaviour of the unit is noted, discontinue use immediately and refer to a qualified EMC engineer.

Contact Address

Digitimer Limited

37 Hydeway
Welwyn Garden City
AL7 3BE
UK

Telephone:- (UK) 01707 328347 (Int.) +44 1707 328347

Fax:- (UK) 01707 373153 (Int.) +44 1707 373153

E-mail:- sales@digitimer.com or technical@digitimer.com

Website:- www.digitimer.com


Servicing & Maintenance

This equipment does not require any regular maintenance but if you would like your DS4 to be serviced we are happy to do so. Please contact us for a reference number and instructions before despatching the unit.

Before each use - The case and all interconnecting cables should be inspected for any damage. The equipment (or the lead) should be sent for repair if any damage is found.

Environmental Considerations

The European Union has adopted Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE), with requirements that went into effect August 13, 2005. WEEE is intended to reduce the disposal of waste from electrical and electronic equipment by establishing guidelines for prevention, reuse, recycling and recovery.

As part of our legal obligation, Digitimer Limited is a registered EEE producer. Our WEEE registration number is WEE/BJ0052TQ. For further information relating to the correct method of disposal of any of our equipment, which features this label , please contact us.

Product Registration

For a speedy response to all your questions now and in the future, please take time out to register your new DS4 at www.digitimer.com/register now! Product registration permits us to quickly advise you of any safety matters or new product information. This web address is your point of contact for all questions regarding the DS4. The site's contents are now growing rapidly, so please bookmark it so that you visit it regularly to check out the new items.

Why Register your purchase?

- Digitimer Limited periodically offers enhancements and updates to our products. Without product registration, users of our products may miss out on announcements of important enhancements to the products that they are using
- Digitimer Limited does not make our customer list available to anyone else. Any information that you provide to us is strictly confidential.

How to Register your purchase

Product registration can be accomplished in two ways. You may fill out and mail in the product registration/warranty card supplied with each Digitimer Limited product. You may also register on-line at our www.digitimer.com/register website.

Product announcement mailing list

Digitimer Limited has e-mailing lists which we use as our primary outlet for announcements of new products, product enhancements and product updates. We strongly recommend that all users of our products sign up for the list that is most appropriate to their area of interest. E-mail is kept to a minimum and list membership is kept in the strictest confidence. Only Digitimer Limited can send mail to members of our e-mailing lists.

You may join the DS4 mailing list through our www.digitimer.com/register website.

Unpacking

After unpacking the DS4 Bi-phasic Stimulus Isolator from the shipping carton, please inspect it for any sign of shipping damage. Please contact the carrier and your Distributor, or Digitimer Limited, immediately if there is any damage. Do not dispose of the shipping carton. The carrier will want to examine the shipping carton to process a damage claim. Digitimer Limited and their Distributors insure all shipments to cover shipping damage.

It is also advisable to keep the shipping carton in the event that the instrument needs to be returned for service.

Included Items

- DS4 Bi-phasic Stimulator
- $\pm 15V$ DC Power Supply (input requirements 110-220VAC)
- Mains lead
- Pair of output plugs (NL985P)
- Operator's Manual (this document)

Optional Accessories

Digitimer can supply extra sets of output plugs (NL985P) for user connection to stimulating electrodes, as well as 1.5m BNC cables (D185-TC3). We also stock sets of the ten batteries (DS4-BATT) which provide the

output power source for the stimulator. One or two stimulators may be mounted in a 19" rack using a specially fabricated frame (model D121-11).

Declaration of Conformity

Ref: N:\Docs\Compend\02\14\Regulatory\CE MARKING\DS4 DoC.doc

DECLARATION OF CONFORMITY

according to ISO/IEC Guide 22 EN45014

Manufacturer's Name: Digitimer Ltd

Manufacturer's Address: 37 Hlydeway
Welwyn Garden City
Hertfordshire,
AL7 3BE
England

as the manufacturer of the apparatus listed, declare that the product (s):

Product Name: Biphasic Constant Current Stimulus Isolator
Model Number: DS4
Product Options: All

conforms to the following Standards:

BS EN 61000-6-1:2007;
BS EN 61000-6-3:2007;
BS EN 61010-1:2001.

The products herewith comply with the requirements of:

the EMC directive 2004/108/EEC
the LV directive 2006/95/EEC

The products were tested in a typical configuration.

Issued on: *March 19, 2009*

Authorised by: 

John R Smale
Technical Director

Functions & Features

Overview of the Controls/Connections

Power Connection

The DS4 control circuitry is powered by the supplied $\pm 15V$ DC power supply adaptor which must be connected to the DS4 via the socket illustrated in Figure 2. The supplied mains lead should then be used to connect the DC power supply adaptor to a suitable mains supply (110-220V AC). Once all connections have been made, the mains socket and the DS4 can be switched on.



Figure 2 Positioning of the DC power supply socket on the DS4.

Note that the DS4 is NOT a re-chargeable unit and the DC power supply must remain connected and switched on while the DS4 is in use.

Power Zone

The DS4 has a POWER switch in the middle of the front panel which switches the unit on or off (by disconnecting the DC power supply). When the power is switched on (toggle switch up), the green LED is lit. It is recommended that the DS4 is switched off when not in use, as significant battery drain will occur if a voltage continues to be applied at the INPUT, when the DS4 is switched on.



Figure 3 The POWER zone of the DS4 front panel.

Gate Zone

A BNC socket is provided to allow a TTL command to "Gate" the output on or off. This feature means that a command voltage signal can be applied continuously at the INPUT and the stimulus turned on and off via digital control at the GATE input. Further, a single command voltage can be applied to multiple DS4 units with individual digital lines used to turn specific DS4 outputs on or off. A toggle switch is provided which enables (UP) or disables (DOWN, labelled OFF) the stimulus output. The amber LED is lit when both the GATE is switched on (toggle switch up) and a TTL high input is being applied at the GATE input (or this input is open circuit), indicating that the output is enabled. This LED will extinguish when the stimulator output is disabled by the Inactivity Sensor. The maximum voltage that can be applied at the GATE input is $\pm 15V$.




Figure 4 The GATE zone of the DS4 front panel.

Input Zone

A command voltage waveform is applied at this BNC socket and is used by the DS4 to determine the shape of the stimulus it delivers. The DS4 accepts various input voltage ranges, but is factory set to $\pm 10V$. An internal jumper allows the operator to change the input voltage range to $\pm 1V$, $\pm 2.5V$ or $\pm 5V$ to match the output range available from the signal generator being used. Please refer to [Internal Jumper Selection](#) for further details. If the input voltage range is changed, this should be indicated by placing a marking in the appropriate box on the front panel. The maximum voltage that can be applied at this input is $\pm 12V$.

The INPUT zone of the DS4 front panel includes a pair of amber LEDs which illuminate when the input voltage exceeds $+0.15\%$ or -0.15% of the full scale voltage.




Figure 5 The INPUT zone of the DS4 front panel.

Output Zone

The OUTPUT zone of the DS4 front panel includes a pair of 2mm touch proof sockets and a range selector switch which allows the operator to adjust the output current range from $\pm 10\mu A$ to $\pm 10mA$, from the compliance voltage source of $\pm 44V$. If the DS4 is set to an input voltage range of $\pm 10V$ and the output range dial is set to $\pm 10mA$, then a $+10V$ signal at the input will deliver $+10mA$ at the output and $+5V$ will deliver $+5mA$ etc.



Figure 6 The OUTPUT zone of the DS4 front panel.

Stimulation electrodes should be connected to the DS4 via the pair of NL985P plugs which are supplied with the DS4.

Output Characteristics

As with use of any electrical stimulator, the operator should be familiar with Ohm's Law ($V=IR$) and check that the stimulation pathway is of low enough impedance to allow the required current to pass through the preparation. For example, if a stimulus of up to $\pm 10\text{mA}$ is required, then the compliance voltage $\pm 44\text{V}$ would be able to pass this current through an impedance of no more than 4.4kohms .

Bear in mind that electrodes in tissues do not behave as resistors - they may have much lower effective series resistance for narrow pulses than for DC inputs, due to polarisation (most metal electrodes polarise with DC inputs). A typical tungsten microelectrode, for example, has a DC resistance of about 200Mohms , an 'impedance' of 1Mohm for 1kHz sine input, and an impedance of perhaps only 100kohm for a $100\mu\text{s}$ pulse, a 2000:1 ratio. Figure 7 predicts that this electrode will pass about $400\mu\text{A}$ for short pulses, but only about $0.2\mu\text{A}$ for very long pulses.




Figure 7 Graph of deliverable current (hatched area) depending upon load.

Checking Actual Output Current

If there is some doubt that the load resistance is low enough to allow operation within the shaded portion of Figure 7, the actual output current can be monitored by inserting a 1kohm resistor into the stimulation pathway as illustrated in Figure 8.

Note that the 1kohm monitoring resistor is connected directly to the preparation ground; this obviously reduces the effective isolation of the DS4 output, which is not normally desirable. All of the output current must however pass through both the 1kohm resistor and the electrode (load resistance), giving an accurate measure on the oscilloscope of the actual output current. The oscilloscope input cannot be put directly across the isolator output because it would then shunt the electrode. The 1kohm series resistor does not contribute appreciably to the total load resistance for high impedance electrodes such as microelectrodes; gross electrodes (e.g. EEG electrode, silver hooks, etc.) may require a smaller series resistor (e.g. 10 ohms) when high output currents are used.




Figure 8 Method to monitor current actually delivered.

Output Capacitance and Risetime

The DS4 has a total output capacitance of approximately 40pF. Thus a 1Mohm resistor across the output will limit the output risetime to about 40 μ s. As discussed above, however, metal microelectrodes which are nominally characterised as having an impedance of 1Mohm (at 1kHz) behave more like 100kohms resistors for short pulses and therefore usually produce risetimes of the order of 10 to 15 μ s. Non-polarisable electrodes (chlorided silver electrodes, micropipettes) behave more like resistors than ordinary metal microelectrodes do; i.e. their impedance is more constant as a function of frequency. Thus, a 10Mohms micropipette acts like a 10Mohms resistor for short pulses as well as long ones. The risetime may therefore be quite long for these electrodes (approximately = 40 μ s electrode DC resistance in Mohms).

How Constant is the Constant Current Output?

The DS4 has more than 900Mohms output resistance at its output terminals. Thus, load resistance variations from a short circuit (0 ohms) to 100Mohms will result in only approximately a 10% variation in output current (assuming that the voltage drop through the load does not exceed the 44V limit discussed above). Obviously few load situations will vary to this extent, and for all practical purposes, the output is constant.

Frequency Response

Figure 9 illustrates the frequency response characteristics of the DS4 when a $\pm 5V$ input is applied with the DS4 set to a maximum input range of $\pm 10V$, an output range of $\pm 10mA$ and a load of $1k\Omega$. The graph shows that the $-3dB$ roll-off is at $50kHz$, while the expected DS4 output is constant for frequencies up to $5kHz$.




Figure 9 Frequency Response Characteristics.

Inactivity Sensor

One of the problems with stimulators that make use of an external voltage source to define a stimulus waveform is that small offsets or noisy baseline signals from the DAC's used to drive them can result in unwanted battery drain or perhaps worse, low amplitude stimulation. The DS4 uses a special "inactivity sensor" to monitor the input voltage and disable the DS4 output if this voltage falls below 0.15% of the full scale value (in -ve or +ve directions) for a user selectable timeout of 200ms, 500ms, 1s or 2s. Unlike other devices which only ever produce an output when the input voltage exceeds a threshold value, this "inactivity sensor" reduces battery usage and damaging "leak currents" during infrequent stimulation, while at the same time maintaining low levels of zero crossing distortion for repetitive waveforms.

The inactivity sensor timeout can be adjusted between the four possible settings using on board jumpers which are accessible from the battery compartment. Please refer to [Internal Jumper Selection](#) for further details.

It is important to note that the inactivity sensor will result in a slight "glitch" in the stimulus waveform when the DS4 is first used after a period of inactivation. This is because the inactivity sensor prevents any output until the command input voltage exceeds $\pm 0.15\%$ of the full scale.




Figure 10 Effect of the inactivity sensor on stimulus initiation.

The example illustrated in Figure 10 shows how the first upward section of a saw tooth waveform (or any similar slowly changing voltage) will have a small distortion at the initiation of the waveform (indicated by the purple circle), but that this distortion is absent in the second and subsequent cycles because the inactivity sensor has not been re-enabled. Clearly, if the slope of the input voltage waveform is very shallow and the timeout setting for the inactivity sensor is short, it is quite possible that the inactivity sensor could reactivate and the “glitch” could repeat. This may be prevented by increasing the duration of the timeout setting or using a steeper waveform shape.

Internal Jumper Selection

The DS4 has two sets of jumpers for (i) input voltage range and (ii) inactivity sensor timeout that can be set by the operator. To gain access to the battery compartment, it is first of all necessary to remove the two screws which hold the cover in place and lift it away from the body of the stimulator. On the right hand side of the battery compartment are two sets of jumpers. The upper set controls the inactivity sensor timeout, while the lower set defines the input voltage range. The battery compartment cover (see Figure 12) includes labelling which indicates jumper positioning.




Figure 11 Location of the internal jumpers.

Battery Life

The ten batteries in the DS4 discharge ONLY when a voltage is applied at the input and battery drain between commanded stimuli is negligible due to the presence of the Inactivity Sensor.

These ten batteries are divided into two sets of five for each polarity phase (a single "Control" battery and four "Stimulus" batteries). The life of "Control" batteries is determined solely by the amplitude and duration of the input pulses and not the output amplitude: continuous operation with high duty cycles at the maximum input amplitude, produces greater battery drain than low duty cycles at low input amplitudes. The drain on this battery is the same for each output range.

Input Voltage (V)	Duty Cycle	Life (hours)
10	50.00%	25
10	10.00%	125
10	1.00%	1250
5	50.00%	50
5	10.00%	250
5	1.00%	2500
1	50.00%	250
1	10.00%	1250
1	1.00%	12500

Table 1 Expected "Control" battery life.

The four "Stimulus" batteries simultaneously discharge at a rate that depends only on the magnitudes of the output pulses. The table below gives a rough guide to the service life for these batteries.

Pulse Amplitude	Duty Cycle	Life to 80% (~34.4V output)	Life to 50% (~22V output)
10mA	50.00%	8 hours	16 hours
10mA	10.00%	40 hours	80 hours
10mA	1.00%	400 hours	800 hours
1mA	50.00%	80 hours	160 hours
1mA	10.00%	400 hours	800 hours
1mA	1.00%	4000 hours	8000 hours

Table 2 Expected "Stimulus" battery life.

As the "Stimulus" batteries discharge, their output voltage decreases. The first number in the Life entries above indicate the service at which the battery voltage in each battery has fallen to 80% (giving a total output voltage available of about 34.4V); the second number indicates the hours of service to battery

voltages of about 50% (total output voltage available of about 22V).

Thus, with fresh batteries, the DS4 will be able to deliver ±10mA through any load resistance less than 4.4kohms, but after 8 hours of continuous service with a 50% duty cycle at the maximum amplitude setting on the 10mA output range, only 7.8mA can be delivered through the 4.4kohms load (i.e. $34.4V/4,400=0.0078A$). In situations where high currents are required and load resistances are near the limits, careful attention must be paid to battery checks and replacement.

In most applications, the entire set has a life of several months, perhaps even approaching the battery shelf life. As the cost of these batteries is relatively low, it is wise to replace all 10 at the same time.

Battery Testing

The DS4 is equipped with battery test sockets which allow the operator to test the status of the batteries without removing them from the unit. Battery voltages can be tested using a digital multimeter, by placing the probes into the labelled test sockets for the +ve and -ve phase and control batteries.

It is important to test the batteries while the DS4 is delivering a long stimulus, so that failing batteries can be more readily identified.

Battery Replacement

If replacement is necessary, the ten batteries are located under the cover panel in base of the unit. Removal of the two screws in the panel will give access.

It is imperative that the batteries are fitted with the correct polarity orientation. This is marked on the battery holders and the cover panel (see Figure 12). The symbols must match.

It is recommended that all ten batteries are replaced at the same time. Sets of replacement batteries are available from Digitimer, using the part code **DS4-BATT**.



Figure 12 DS4 battery compartment labelling.

Rack Mounting

The DS4 can be mounted and positioned on any plane and 19" rack mounting is possible by using two of the box fixing screws (diagonal corners) through our D121-11 Mounting Frame, which is available from Digitimer or our representatives. The mounting frame allows any two of our DG2A, DS2A, DS3 or DS4 devices to be positioned side by side.

References

As the DS4 is a brand new product we do not yet have any publications that cite its use. However, if you publish research which has used the DS4, please cite it in your methods section to help other researchers. Digitimer would appreciate a copy of any relevant publications and could add details to this section of the manual.

DS4 Specifications

Output:	Bi-phasic constant current proportional to the input voltage
Output Ranges:	$\pm 10\mu\text{A}$; $\pm 100\mu\text{A}$; $\pm 1\text{mA}$; $\pm 10\text{mA}$ for a full scale input
Output Duration:	$>2\mu\text{s}$
Compliance:	$\pm 44\text{V}$ from 8x GP23A batteries
Linearity:	$\pm 3\%$ of full scale output for each output range
Output Impedance:	$>900\text{Mohms}$
Output Rise Time:	$<5\mu\text{s}$ (1kohm load), $<40\mu\text{s}$ (1Mohm load)
Inputs:	
IN:	Ranges: ± 1 ; ± 2.5 ; ± 5 ; $\pm 10\text{V}$ full scale (selected by an internal jumper) with a limit of $\pm 12\text{V}$ max. without damage. Input Impedance: 1Mohm
GATE:	Range: TTL; Gate OFF if Low; Gate ON if High or open circuit. Limit of $\pm 15\text{V}$ max. Input Impedance: 10kohm
Inactivity Sensor:	The output is disabled if the voltage input remains below 0.15% of the full scale value (-ve or +ve phases) for a user selectable period of 200ms, 500ms, 1s or 2s. This time period can be adjusted with an internal jumper.
Connections:	Output - 2mm shrouded, touch-proof sockets (red and black) spaced at 0.75" Input - Front panel BNC socket Gate - Front panel BNC socket Battery Test - Six 2mm sockets Power - Socket for external power supply
Controls:	Gate - On/Off toggle (Off overrides BNC input) Output Range - 4 position rotary switch Power - On/Off toggle switch
Indicators:	
Power ON	LED Green (lit when the power supply is connected and DS4 is switched On)
Gate Enabled	LED Amber (lit when Gate is On and the Gate Input is held TTL high or open circuit)
Phase +ve	LED Amber (lit when input exceeds +0.15% of full scale voltage)
Phase -ve	LED Amber (lit when input exceeds -0.15% of full scale voltage)
Power:	Included external power supply (input voltage 110V - 220V) providing $\pm 15\text{V}$ DC output. 10 x 12V GP23A Batteries.
Mounting:	One or two stimulators may be mounted in a 19" rack using a specially fabricated frame (model D121-11) available from Digitimer Ltd.
Dimensions:	190 x 110 x 80 (w x h x d)
Weight:	500g (approx.)

Warranty

Digitimer Limited warrants to the first purchaser, for a period of one year from the date of purchase, that this Digitimer Instrument (hereafter referred to as the "Product") will be free from defective workmanship and materials, and agrees that it will, at its option, either repair the defect or replace the defective Product or part thereof at no charge to the purchaser for parts and labour. The Product must be returned to Digitimer Limited, carriage paid and insured.

Digitimer Limited will return the Product, carriage paid and insured, in the most appropriate method as determined by Digitimer Limited. If a faster shipping service is desired by the customer, any additional special delivery expenses must be paid by the customer.

This warranty does not apply to shipping damage. Digitimer Limited fully insures all shipments. Any claims of damage upon receipt must be filed with the carrier and Digitimer Limited immediately. This warranty does not apply to any exterior appearance item of the Product which has been damaged or defaced, which has been subjected to misuse and abuse, abnormal service or handling, or which has been altered or modified in design or construction.

This warranty does not apply to any interconnection cables supplied with the Product.

This warranty does not apply if any unauthorised repairs, modifications or alterations have been made to the Product.

This warranty applies to software products only to the extent of maintenance release software to correct improper operation of the Product. Software updates to increase the capabilities of the present product are not to be provided under the terms of this warranty. Updates will be sent at no cost to the customer by normal common carrier routes. If faster delivery is desired, the customer must assume any additional expenses for special delivery service.

No sales organisations, other than Digitimer Limited itself, are authorised to make any warranties other than those described above, or to extend the duration of any warranties beyond the time period described above on the behalf of Digitimer Limited. If Digitimer Limited agrees to such a modification of this warranty, Digitimer will furnish a modified copy of this agreement, which must be presented if a claim is being made under these modified terms.

Obtaining Warranty Service

Warranty service of this Product can be obtained by returning the Product, carriage paid and insured, to Digitimer Limited, or the Distributor from whom the equipment was purchased. Prior authorisation before shipping the product is advised for the most expedient service.

Product change or discontinuation

Digitimer reserve the right to discontinue any instrument or to change its specification without notice, and without responsibility for incorporating changes in instruments already sold.

Operator Notes

Please use the space below to make notes on particular settings or include any warranty information or service history for your DS4.

Digitimer Limited

37 Hydeway
Welwyn Garden City
AL7 3BE
UK

Telephone:- (UK) 01707 328347 (Int.) +44 1707 328347

Fax:- (UK) 01707 373153 (Int.) +44 1707 373153

E-mail:- sales@digitimer.com or technical@digitimer.com

Website:- www.digitimer.com