

"The EasyCode software you provided is simple to use. We use a large screen Macintosh for creating and storing protocols. We especially like the ability to create new protocols and download them directly to the ValveBank, it really couldn't be easier."

**Dr. Susan Abrahamson**  
University of California,  
Berkeley


### EasyCode Dialog Box

- Enter start and stop times accurate to 0.01 sec.
- Use all ValveBank commands including multi-channel looping and interface triggering

### Free your computer during experiments

Transfer and run programs stored in the ValveBank, which can be disconnected from the computer – leaving it available for data acquisition.

# Expand the computing power of your ValveBank system.



### EasyCode® Software

Program your ValveBank® with a Macintosh or PC using easy "click-and-drag" time bars. Save and load unlimited programs to disk, print out program listings, then download your sequences into the ValveBank's memory with the included serial cable in one simple step.

- **Program up to sixteen channels of valves and digital outputs**  
Open multiple experiment windows. List sequences to screen or printer. Copy and paste. Zoom in and out of your experiment.

### Ordering Information

Part No.	EasyCode Software
01-02-USB	EasyCode® Mac - ValveBank programming software Macintosh OS 6-9 & Classic with USB-to-serial adapter
01-03	EasyCode® PC - ValveBank programming software - Windows 3.1-95
01-03-2000	EasyCode® PC - ValveBank programming software - Windows 98-XP
USB	USB to serial converter - PC or Macintosh

U.S./Canada prices shown. International prices add 15%. Email or visit web store for latest prices.

