

# Industry-standard electrical noise protection.



## Type II Faraday Cages

The Type II Faraday Cage offers improved access and simplified assembly. The “window-shade” type retracting front panel is easier to operate than hinged doors and causes less disturbance when adjusted. This front panel may be positioned anywhere between fully opened and closed and stays in position without a fastener. The front door is shipped assembled and the entire unit may be assembled in a few minutes with a screwdriver (provided).

This cage incorporates the same stainless steel frame and copper-mesh material as previous versions. It mounts to (and requires) TMC’s full-perimeter enclosures and mounts to our 63-500 Series tables.

## Ordering Information

Part No.	TMC Air Table
TM-81-333-03	Type II Faraday Cage, 40 in.
TM-81-334-03	BenchTop Faraday Cage
	Complete TMC product line available. Call or email for part numbers.

U.S./Canada prices shown. International prices add 15%. Email or visit web store for latest prices.

800.998.MATE | www.autom8.com | 812 Page Street, Berkeley, CA 94710 USA  
tel 510.845.6283 | fax 510.280.3795 | e-mail info@autom8.com

## BenchTop Faraday Cage

Until recently, TMC's Faraday Cages have only been available as part of a complete vibration isolation table system. The cages required our Perimeter Enclosure option and could not be assembled without a TMC 63-500 Series Vibration Isolation Table as its base.

We now offer the same line of 40 in. tall cages with a baseplate which allows the cage to be used on a bench-top without a corresponding TMC table. The base of the cage is a reinforced stainless steel plate which can support a compact vibration isolation system, microscope, or other instrument.

## Options and Features

### Type II Cages

- 2 in. diameter holes for cable passage
- Optional armrest pads are now available with Faraday Cages
- Optional "U" shaped hanging shelf
- Optional non-isolated sliding shelf oriented front-to-rear
- Full perimeter enclosure required to mount cage

### BenchTop Faraday Cage

- BenchTop cage includes stainless steel baseplate
- Bench-mounted cage does not require vibration table

