

Integrated Digitizer

- 40 kHz sampling rate
- 16-bit A/D (18-bit internal resolution)
- Stimulus voltage ranges: ± 250 mV and ± 2000 mV
- Zap voltage range of ± 1000 mV



Feedback Gain Settings

- 10 M Ω , 100 M Ω , 1 G Ω , 3.3 G Ω , 10 G Ω

Low RMS Noise (DC to 3kHz)


- 0.3pA @ 10G Ω , 1.0pA @ 1G Ω , 7pA @ 100M Ω

Compensations

- Up to 4 Cap. Compensations
Cfast x 1, Cslow x 3
- 0-100 pF per compensation
- Series Resistance Compensation
Offset Compensation (± 250 mV) Optional Active Leak Compensation

Current Clamp

- ± 2 nA range with 1.25 pA res.
- ± 20 nA range with 12.5 pA res.
- ± 200 nA range with 125 pA res.



The patented Spread Frequency Compensation can automatically compensate any arbitrary capacitance profile in approximately 3 seconds.

The smallest patch clamp system available...



USB-powered Miniature Patch Clamp Amplifier and Digitizer

Full-featured 1-channel amplifier

Pico is a feature-rich, low-noise patch clamp amplifier with an integrated digitizer and headstage. Pico is ideally suited for whole cell (V_{clamp} and I_{clamp}) and single channel (patch, planar lipid bilayer and synthetic nanopore) recording, as well as cellular electrochemistry.

Voltage Clamp and Current Clamp

Clean Head Switching™ technology allows software-controlled switching between voltage clamp and current clamp, or between multiple feedback resistors, without introducing any artifacts. This enables one Pico to support a wide range of applications including whole cell, single channel, multi-cell, bilayers, electrochemistry, and current clamp.

Software

Compatible with: WinWCP by the University of Strathclyde. jClamp by SciSoft Company. TecellaLab software, with Data Export to ATF, tab-formats. SDK/API available. Other 3rd party software support.

Ordering Information

Part No.	Pico 2
Te-Pico 2	Pico 2 Mini electrophys. amp, BOB, and digitizer
Te-Trigger Cable	Optional cable to control external devices with WinWCP
Te-jClamp32	jClamp Pico Software academic license (call for non-acad)



AUTOMATE SCIENTIFIC